


PORTABLE SOLUTIONS for Gas, Liquid & Solid Determination with Lab Quality Confirmation

Excellent for All First Responders, Forensics & Scientists

Limited Time Offer!

VIR - Gas Analysis Detection System

For a limited time, JASCO is offering the VIR GAS System with a complimentary Dell Latitude ATG rugged laptop & pelican case for safe transport & storage



The VIR gas analysis detection system checks gases detected by a fourier transform infrared spectrometer against a library of over 5000 gases. This system has the ability to be automated from inside a vehicle, limiting first responders exposure to outside air prior to analysis.

Common Uses:

- First Response Vehicles
- Firefighting & Rescue Vehicles
- Monitoring Services



Features Include:

- Automated system
- Can accurately check air every 15-30 seconds
- New Gas spectra can be added by users for proprietary data or unique gas compounds not currently in the database


Monitoring & Automated Reporting:

- Can be triggered to monitor gases in assigned intervals in many environments
- Fast start-up
- 120 calibrations included (more optional)
- Automated report includes detected gases and quantitative amount (if calibration exists)

VIR with IRT-1000 Microscope

For a limited time, JASCO is offering the VIR Microscopy System with a complimentary Georgia State Crime database & pelican case for safe transport & storage

The integration of the Irtromu with the VIR-100 portable FT-IR spectrometer provides a portable FT-IR microscopy system compact and light-weight for field operation. This system enables FT-IR microscopy analyses to go from daily use in the lab, directly into the field.



System Benefits:

- Systems offer spectral resolution and sensitivity needed for accurate QC, Identification & Quantitation
- Provides a dependable, accurate solution
- Sample examination at site of incidence, minimize risk of contamination
- Save time and money. Portable solution that can be used at multiple facilities
- Minimal training, simple operation

Features Include:

- Fast sample data acquisition
- Quick instrument set-up
- Validation software included
- Optional Micro ATR is available for Irtromu
- Optional 2nd detector position for MCT detector
- Can also be used with ATR, diffuse reflection, or transmission accessories

